

Product: Egg Yolk Tellurite Emulsion 20%

Catalog Number: iS47021

Lot Number:

Date of Release:

Expiration Date:

Product Description: Sterile stabilized tellurite emulsion of egg yolk recommended for identification of the Staphylococcus species.

Components:

Egg Yolk	30.0 mls
Sterile Saline	64.0 mls
Sterile 3.5% Potassium Tellurite Solution	6.0 mls

Preparation Instructions: Warm the refrigerated Egg Yolk Tellurite Emulsion to 45-50 °C and aseptically add 50 ml to 950 ml of sterile, molten Baird Parker Agar Base. Mix well and pour into sterile petri plates.

Storage Temperature: Store at 2-8 °C.

Precautions and Disclaimer: For laboratory use only. Not for drug, household or other uses.

QC Department